

# VariLoc

Variable Angle Locking Technology




THERAPY OPTIMIZED  
**EXPECT EXTRAORDINARY**

## Variable Angle Locking Technology


### Mechanism

- The thread of the Variable Angle screw miss-match the thread of the locking hole by 1 pitch to achieve certain degree of off-axis fixation.




### Advantages

- Unique variable angle locking screws enable the surgeons to adapt the screw angulations to patient anatomy



10°-2.4/2.7/3.5/5.0mm VA screw

- Special Instrument Facilitates the Insertion of VA Screw



- Diverse Locking Screw Options

